


# Tips to Improve NGS Library Preparation



## Tips to Improve

### **NGS Library Preparation**

A high-quality NGS library is essential in the generation of reliable next-generation sequencing data. Efforts in improving your NGS library construction pipeline is time well spent, as following established best practices and optimization of key elements along the process will yield reproducible results that you can count on. Depending on the goals of your lab, various kits are available for different types of analysis and some are more amenable to high-throughput automation. Navigating potential pitfalls that necessitate "re-dos" can reduce costs, save time, and reduce overall frustration in the lab. These tips to improve your NGS library prep will help you to get to high-quality libraries and better data in no time.




#### High quality DNA/RNA In

To create an NGS library it is critical to begin with a pure and highly concentrated  $\ensuremath{\mathsf{DNA}}$  or RNA sample. Make sure to try various extraction approaches and choose the solution that yields the purest DNA or RNA with sufficient input for several runs. Further, choose a library preparation kit that has been developed to support the target analysis, starting material, and the degree of automation that you intend to establish in your laboratory.



#### Concentration and cycle optimization of PCR adapters

During the library preparation process, target fragments are amplified to allow for sufficient material for the degree of coverage required to generate a trustworthy sequence determination. However, over-amplification can lead to PCR bias, duplicates, and adapter dimers. Optimizing the number of cycles based on the input concentration can reduce these unintended biproducts leading to a more efficient and higher quality analysis.



#### Library QC and quantification

The approach to NGS library QC is well established and includes both size determination and quantification of libraries. Accurate and consistent sizing and concentration of libraries indicates high-quality library preparation and usually guarantees reliable coverage and reproducible data in sequencing experiments.

Sizing via microfluidic or capillary electrophoresis will show libraries that are not of high enough quality which will allow you to react early rather than wasting precious capacity on the flow cell with fragments that are not of sufficient quality for sequencing. Quantification using an intercalating dye and fluorescence detector will allow for library normalization so that the same mass of each sample is evenly distributed and represented in the data set. Monitoring this data over time can uncover variability in your library construction pipeline and provide key insights that can be used to optimize your manual or automated method for preparing libraries.



## 4.

## Automation and the Power of Standardized and Tested Protocols




Although most liquid handling platforms can be programmed to generate libraries, many vendors treat library preparation as a custom project requiring solution configuration and weeks of software programming and testing. Choose an automation provider that offers a standardized configuration for NGS and a library of standardized and tested protocols that can be available within days of installation. There is tremendous power in using protocols that are consistently utilized by researchers from around the globe on the exact hardware configuration that you are investing in.

For instance, every customer that purchases a G.STATION from DISPEN-DIX receives exactly the same hardware configuration and can choose from a library of protocols in the cloud that have been developed in house by experts with decades of experience optimizing and troubleshooting NGS Library Preparation pipelines. This will save you months of effort in optimizing and developing a method in house, and provides a welcomed piece of mind when making such a substantial investment in automation.

## 5.

#### Optimize for Costs where it Makes Sense



Although the cost of high-quality library preparation kits and quality control are always recommended, there are some aspects of the library preparation process that do not need to be expensive. Automation consumables, especially pipette tips, can be a cost burden to the lab and can be scarcely available from vendors. Look into new technologies that enable tip-free dispensing and bead-based clean-ups to optimize for cost per sample.

The I.DOT non-contact dispenser reliably dispenses from 8 nL to 10s of uL enabling tip-free enzyme, buffer, and bead addition as part of a semi-automated or fully automated workflow. The C.WASH Bead-Based Clean-up Device allows for ethanol addition and removal in minutes with zero consumable usage. These technologies have been designed for NGS and offer a much more pragmatic solution than traditional tip-based liquid handlers.







©2021 BICO AB. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of BICO is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO may refer to the products or services offered by other companies by their brand name or company name solely for clarity and does not claim any rights to those third-party marks or names. BICO products may be covered by one or more patents. The use of products described herein is subject to BICO's terms and conditions of sale and such other terms that have been agreed to in writing between BICO and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of BICO products in practicing the methods set forth herein has not been validated by BICO, and such nonvolidated use is NOT COVERED BY BICO'S STANDARD WARRANTY, AND BICO HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner BICO's terms and conditions of sale for the instruments, consumables or software mentioned, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by BICO that it currently or will at any time in the future offer or in any way support any application set forth herein.

#### Contact

Tel: +49 (0) 711 490 544 00 Email: info@dispendix.com Website: www.dispendix.com